

-53-

again Minipreps for Sequencing

MINIPREP FOR DNA PREPARATION

Wd: 1, 3, 4, 5, 6, 10,

NCW: 13, 15, 16, 17, 18, 19, 20

- 1) INOCULATE CLONES INTO 3.0 ML OF LB MEDIA OVERNIGHT AT 37 DEG.
 - 2) THE FOLLOWING DAY LABEL 15 ML. STERILIZED EPPENDORF TUBES TO BE USED. *2ml*
 - 3) FILL 3/4 OF EACH TUBE WITH THE CULTURED BACTERIA. CLOSE THE CAP
 - 4) SPIN TUBES USING THE EPPENDORF CENTRIFUGE FOR ONE MINUTE AT 6000 RPM.
PREPARE THE POT FOR BOILING WATER
 - 5) DECANT THE SUPERNATANT, LEAVING A LITTLE BIT IN THE TUBE.
 - 6) RESUSPEND THE PELLET BY VORTEXING
 - 7) ADD 120 micro ml OF STEL-T SOLUTION. MIX GENTLY BY INVERTING THE TUBES. DO NOT VORTEX. *240ul*
 - 8) KEEP THE TUBES IN BOILING WATER FOR EXACTLY 1 MINUTE AND THEN IMMEDIATELY TRANSFER THE TUBES TO A BUCKET OF ICE.
After boiling the DNA appears whitish in colour.
 - 9) KEEP THE TUBES IN ICE FOR 5 MINUTES. *→ ic* *61 + 2ml RNase (1.10.6-2) [5-10]*
 - 10) SPIN THE TUBES IN THE MICROFUGE FOR 5 MINUTES AT THE HIGHEST SPEED. *(12 x 1000 R.M.)*
 - 11) USING A STERILIZED TOOTHPICK PICK UP THE PRECIPITATE AND DISCARD IT, LEAVING THE SUPERNATANT BEHIND.
 - 12) ADD 120 micro ml OF 2-PROPRANOL (same quantity as the STEL-T solution) 2-ISOPRANOL PRECIPITATES THE DNA. MIX THE TUBES BY INVERSION AND KEEP ON ICE FOR 5 MINUTES. *240ul*
 - 13) SPIN FOR 5 MINUTES IN THE MICROFUGE AT THE HIGHEST SPEED. *(12 x 1000 R.M.)*
 - 14) DECANT THE SUPERNATANT AND LEAVE THE TUBES TO AIR DRY IN AN INVERTED POSITION.
(To hasten this step, you may add 200 micro ml of 100% alcohol without disturbing the pellet and allow to air dry as before.)
discard the alcohol.
 - 15) ONCE THE PELLET IS DRY, WHICH CAN BE MADE OUT BY THE POWDERY APPEARANCE OF THE PELLET, RESUSPEND IN 60 micro ml OF TE WITH RNAase. *170ul*
- WORKING SOLUTION OF TE WITH RNAase
- 2 micro ml of RNAase (20 mg/ml) in 1 ml of TE
- 16) THESE DNA CAN BE FROZEN AT -20 TEMPERATURE IF REQUIRED.

Digest: 25 µl DNA

32 µl 10x Buffer 3

3 µl 10x KCl

0.5 µl NotI

0.5 µl SNI

x14

44.8 µl

42 µl

7 µl
7 µl

Take 7 µl 10 DNA

+ 1 µl RSC + 24 µl H₂O

24 37°C

Ret 100µl → ad 100µl TC

→ 200µl PCIA → 6

→ 200µl CIN

→ 200µl NucAc

→ 50µl Glat 1h -20°C


→ 18 70% Glat

→ ⑦ is gone



sequencing

Not-BU1 digest:



⑧ sequencing 3, 4, 10, 13, 16

-55-

Sequencing:

5' primer Gwars TBOUHO → Schwachman

3' primer CT-BM-SFG primer

Glas's 3, 4, 10, 13, 16

but: USB → Archival "Pentamers"

Respending in 2 μ l H₂O.1 μ l DNA2 μ l 5 \times Sequencase.1 μ l SFG primer.6 μ l H₂O

pe Reaction

3' end

↓
100°C 3-4 min↓
dry ice (15') perhaps looking↓
15 sec spin↓
3ccre Reaction: 1 μ l S-ATP (7000 U/ml)1 μ l DTT 1mM (100 U/ml)2 μ l labeling mix preheated at 5'2 μ l Sequense 1.8 U/ml↓
4 min RT1 μ l H₂O10 μ l70 μ l2 μ l6 μ l

Reaction

3.5 μ l - to 2.5 μ l aliquoted ddH₂O

↓

5' 37°C.

⊕ 4 μ l Stop solution.

→ 5% sequencing Gel:

4.75g	DM
0.25g	Bis
42g	Urea
10mM	10% BGE
500 μ l	APS 10%
40 μ l	Temed

Loading:

5' Reaction 3, 4, 10, 13, 16

3' Reaction 3, 4, 10, 13, 16.

Run: 1322V / 34 μ A / 45 Watts

Exposure 4¹⁰ - 11⁰⁰ → 5h.

Amplifier von Beckman: 1, 3, 4, 5, 6, 10, 13, 16, 17, 18, 19

-52-

Sequences

⇒ 5' reaction equilibrium was, because
TRG-UOH is phasing in CDL-type!

⇒ 3' Reaction:

CCC GGC NAGA / ~~SCG~~ / ~~ATG~~ / ~~GGC~~ / ~~CGC~~ / ~~GAA~~ / ~~GAC~~ / ~~GGT~~ / ~~GCT~~
 GGG ~~CNG~~ TCT / ~~CCG~~ / ~~TAG~~ / ~~CCG~~ / ~~ACA~~ / ~~CTT~~ / ~~CAG~~ / ~~CGA~~ / ~~CCA~~ / ~~CCG~~
 3) ~~GGG~~ / ~~GCG~~ / ~~TAG~~ / ~~GCG~~ / ~~TGG~~ / ~~TAG~~ / ~~TGG~~ / ~~TGG~~ / ~~GGG~~ / ~~GTT~~ / ~~GGG~~
~~CCG~~ / ~~CGA~~ / ~~ATC~~ / ~~CCG~~ / ~~ACC~~ / ~~ATT~~ / ~~ACG~~ / ~~AAZ~~ / ~~GCK~~ / ~~CAT~~ / ~~CCG~~
 CCC GGC NAGA / ~~SCG~~ / ~~ATG~~ / ~~GGC~~ / ~~CGC~~ / ~~GAA~~ / ~~GAC~~ / ~~GGT~~ / ~~GCT~~
 GGG ~~CNG~~ TCT / ~~CCG~~ / ~~TAG~~ / ~~CCG~~ / ~~ACA~~ / ~~CTT~~ / ~~CAG~~ / ~~CGA~~ / ~~CCA~~ / ~~CCG~~
 GTC / ACT / GTC / GAG / GAG / CTT / ~~GAA~~ / ~~CTT~~ / ~~ACT~~ / CTT / CC
 CAG / TGA / CAC / GTC / CTC / ~~GAA~~ / ~~CTT~~ / ~~ACT~~ / ~~ACT~~ / CTT / CC
 GAA ATA TAA GGC NAGA / ~~SCG~~ / ~~ATG~~ / ~~GGC~~ / ~~CGC~~ / ~~GAA~~ / ~~GAC~~ / ~~GGT~~ / ~~GCT~~

4) CCC GTC AGA / ~~SCG~~ / ~~ATG~~ / ~~GGC~~ / ~~CGC~~ / ~~GAA~~ / ~~GAC~~ / ~~GGT~~ / ~~GCT~~
 GGG ~~CNG~~ TCT / ~~CCG~~ / ~~TAG~~ / ~~CCG~~ / ~~ACA~~ / ~~CTT~~ / ~~CAG~~ / ~~CGA~~ / ~~CCA~~ / ~~CCG~~

5) CCC GTC AGA / ~~SCG~~ / ~~ATG~~ / ~~GGC~~ / ~~CGC~~ / ~~GAA~~ / ~~GAC~~ / ~~GGT~~ / ~~GCT~~

6) CCC GTC AGA / ~~SCG~~ / ~~ATG~~ / ~~GGC~~ / ~~CGC~~ / ~~GAA~~ / ~~GAC~~ / ~~GGT~~ / ~~GCT~~

10) CCC GTC AGA / ~~SCG~~ / ~~ATG~~ / ~~GGC~~ / ~~CGC~~ / ~~GAA~~ / ~~GAC~~ / ~~GGT~~ / ~~GCT~~
 GGG ~~CNG~~ TCT / ~~CCG~~ / ~~TAG~~ / ~~CCG~~ / ~~ACA~~ / ~~CTT~~ / ~~CAG~~ / ~~CGA~~ / ~~CCA~~ / ~~CCG~~

11) CCC GTC AGA / ~~SCG~~ / ~~ATG~~ / ~~GGC~~ / ~~CGC~~ / ~~GAA~~ / ~~GAC~~ / ~~GGT~~ / ~~GCT~~

12) CCC GTC AGA / ~~SCG~~ / ~~ATG~~ / ~~GGC~~ / ~~CGC~~ / ~~GAA~~ / ~~GAC~~ / ~~GGT~~ / ~~GCT~~
 GGG ~~CNG~~ TCT / ~~CCG~~ / ~~TAG~~ / ~~CCG~~ / ~~ACA~~ / ~~CTT~~ / ~~CAG~~ / ~~CGA~~ / ~~CCA~~ / ~~CCG~~

13)

16)

Result: all have CDL-type → look 5' Prime!