

-53-

again Minipreps for Sequencing

## MINIPREP FOR DNA PREPARATION

Wd: 1, 3, 4, 5, 6, 10,

NCW: 13, 15, 16, 17, 18, 19, 20

- 1] INOCULATE CLONES INTO 3.0 ML OF LB MEDIA OVERNIGHT AT 37 DEG.
- 2] THE FOLLOWING DAY LABEL 1.5 ML. STERILIZED EPPENDORF TUBES TO BE USED. *2ml*
- 3] FILL 3/4 OF EACH TUBE WITH THE CULTURED BACTERIA. CLOSE THE CAP
- 4] SPIN TUBES USING THE EPPENDORF CENTRIFUGE FOR ONE MINUTE AT 6000 RPM.  
PREPARE THE POT FOR BOILING WATER
- 5] DECANT THE SUPERNATANT, LEAVING A LITTLE BIT IN THE TUBE.
- 6] RESUSPEND THE PELLET BY VORTEXING
- 7] ADD 120 micro ml OF STEL-T SOLUTION. MIX GENTLY BY INVERTING THE TUBES. DO NOT VORTEX. *240ml*
- 8] KEEP THE TUBES IN BOILING WATER FOR EXACTLY 1 MINUTE AND THEN IMMEDIATELY TRANSFER THE TUBES TO A BUCKET OF ICE.  
*After boiling the DNA appears whitish in colour.*
- 9] KEEP THE TUBES IN ICE FOR 5 MINUTES. *→ ic* *61* *+ 2ml RNase* *(1.10.0-2)* *[5.1.10]*
- 10] SPIN THE TUBES IN THE MICROFUGE FOR 5 MINUTES AT THE HIGHEST SPEED. *(12 x 10000 RPM)*
- 11] USING A STERILIZED TOOTHPICK PICK UP THE PRECIPITATE AND DISCARD IT, LEAVING THE SUPERNATANT BEHIND.
- 12] ADD 120 micro ml OF 2-PROPRANOL (same quantity as the STEL-T solution) 2-ISOPRANOL PRECIPITATES THE DNA. MIX THE TUBES BY INVERSION AND KEEP ON ICE FOR 5 MINUTES. *240ml*
- 13] SPIN FOR 5 MINUTES IN THE MICROFUGE AT THE HIGHEST SPEED. *(12 x 1000 RPM)*
- 14] DECANT THE SUPERNATANT AND LEAVE THE TUBES TO AIR DRY IN AN INVERTED POSITION.  
*(To hasten this step, you may add 200 micro ml of 100% alcohol without disturbing the pellet and allow to air dry as before.)*  
*discard the alcohol.*
- 15] ONCE THE PELLET IS DRY, WHICH CAN BE MADE OUT BY THE POWDERY APPEARANCE OF THE PELLET, RESUSPEND IN 60 micro ml OF TE WITH RNAase. *170µl*

## WORKING SOLUTION OF TE WITH RNAase

2 micro ml of RNAase (20 mg/ml) in 1 ml of TE

- 16] THESE DNA CAN BE FROZEN AT -20 TEMPERATURE IF REQUIRED.

Digest: 25 µl DNA

32 µl 10x Buffer 3

3 µl 10x KCl

0.5 µl NotI

0.5 µl SNI

214

44.8 µl

42 µl

7 µl  
7 µl


Take 7 µl 10 DNA

+ 1 µl RSC + 24 µl H<sub>2</sub>O

24 37°C

Ret 100µl → ad 100µl TC  
 → 200µl PCIA → 6  
 → 200µl CIN  
 → 200µl NucAc  
 → 500µl G101 1h -20°C  
 → 18 70% G101  
 ↓  
 sequencing

Not-BU1 digest:



⊗ sequencing 3, 4, 5, 10, 13, 16

-55-

Sequencing:

5' primer Gucens TBOUMO → Schwachman

3' primer CT-BM-SFG primer

Clare's 3, 4, 10, 13, 16

but: USB → Rockwell "Pentamers"

Resuspending in 2 μl H<sub>2</sub>O.

1 μl DNA

2 μl 5' Sequencer

1 μl SFG primer

6 μl H<sub>2</sub>O

pe Reaction

3' end

↓  
100°C 3-4 min↓  
dry ice (15') perhaps looking↓  
15 sec spin↓  
3 sec

re Reaction: 1 μl S-ATP (7000 cpm/μl)

1 μl DTT mix (100 μM)

2 μl labeling mix preheated at 5'

2 μl Sequencer 1.8 labeled

↓  
4 min RT

1 μl

10 μl

70 μl

20 μl

6 μl

Reaction

3.5  $\mu$ l - to 2.5  $\mu$ l aliquoted ddH<sub>2</sub>O

↓

5' 37°C.

⊕ 4  $\mu$ l Stop solution.

→ 5% sequencing Gel:

4.75g	DM
0.25g	Bis
42g	Urea
10mM	10% BCT
500 $\mu$ l	APS 10%
40 $\mu$ l	Temed

Loading:

5' Reaction 3, 4, 10, 13, 16

3' Reaction 3, 4, 10, 13, 16.

Run: 1322V / 34  $\mu$ A / 45 Watts

Exposure 4<sup>10</sup> - 11<sup>00</sup> → 5h.

Amplifier von Beckman: 1, 3, 4, 5, 6, 10, 13, 16, 17, 18, 19

